Tedequim

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TEDEQUIM®

BIOCERAMIC ROOT CANAL SEALER

BASED ON CALCIUM SILICATE TECHNOLOGY

This sealer, whose powder is mainly composed of tricalcium and dicalcium silicates, silicon oxide and zirconium oxide, has excellent biocompatibility, antimicrobial action and bioactivity. The liquid is mainly composed of bidistilled water, calcium chloride and a hydrophilic plasticizing agent (water-soluble polymer).

Due to its physicochemical properties, it continues to harden in the presence of dentin and tissue moisture, as well as blood contamination. It is dimensionally stable, and it expands slightly upon hardening to ensure an adequate seal of the root canal. It remains hard, compact, and insoluble over time.

During preparation, when mixing powder and liquid, its pH is close to 12 due to the release of calcium hydroxide and oxydriyl radicals by the hydration of tricalcium silicate. If during hardening it comes into contact with tissue fluids, the calcium released reacts with the phosphate ions present in the latter and the phosphatase enzymes producing calcium phosphate in the form of hydroxyapatite.

Indication:

Permanent obturation of root canals with gutta-percha cone after endodontic treatment in inflamed or necrotic pulps.
Single cone technique or cold lateral compaction.
Retreatment of root canals.

Technical specifications:

Setting time: between 2 and 3 hours
Working time: 7 to 10 minutes
Fluidity: adequate and according to ISO 6876 for a mixture of 1 measure of powder and 1 drop of liquid.
Radiopacity: between 5 and 6 mm Al
Solubility: 2.8% according to ISO 6876 (2012) requirements.

Characteristics:

- Highly biocompatible
- Bactericidal due to its high alkalinity at the time of hardening, which maintains it over time.
- Good radiopacity.
- Induces bone generation and chemical adhesion to the root dentin.
- Produces a hermetic seal of the root canal due to its slight setting expansion avoiding bacterial filtration.
- Easy to remove from the root canal in case of retreatment.
- Does not color the tooth structure.
- Does not contain resins.
- Low solubility.

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Instructions for use:

- 1- Shake the powder bottle before opening.
- 2- Remove a measure of the powder flush and place it on a glass slab.
- 3- Immediately cap the powder and store it inside the bag together with the moisture absorber of silica gel.
- 4- Shake the liquid well.
- 5- Place a drop of liquid next to the powder.
- 6- Spatulate until a mixture forms threads of 2 cm in height in 45 to 60 seconds.
- 7- Store all the components, well closed in the metallized bag.

Before placing the sealer in the root canal:

- 8- Clean and disinfect the canal according to endodontic protocol. Then rinse with distilled water or physiological solution to remove the remains of the irrigation solutions. Slightly dry the canal.
- 9- Select the appropriate gutta-percha tip.
- 10- Obturate with the single cone technique (recommended) or cold lateral compaction.

Contraindications:

Do not use if the patient is allergic to any of the ingredients.
Do not use in temporary teeth.

Warnings and precautions:

The use of the product is at the discretion of the professional, according to the clinical condition of the patient.
Once opened, the powder should be used within 6 months.
Always shake the liquid before use.

Preservation and storage:

Keep tightly closed, inside the metallized bag together with the moisture absorber at a temperature not exceeding 30 °C.

Composition:

- Powder: white, mainly composed of tricalcium silicate, dicalcium silicate, silicon oxide and zirconium oxide.
- Liquid: Calcium chloride, hydrophilic plasticizing agent and water

Contents:

- Powder: 10 g
- Liquid: 10 ml

Autorized by ANMAT PM 1277-35

Technical director: Miryam Della Vedova - Pharmacist - Licence 2107

Exclusive sale to professionals and sanitary facilities

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